



SDactiv

Surface Disinfectant

CONTAINS ALCOHOL



63%
ALCOHOL
CONTENT

60 SEC
CONTACT
TIME

LEMON
AROMA

STREAK
FREE

KEY FEATURES

Bossklein SDactiv Surface Disinfectant is a powerful, alcohol based disinfectant spray. Designed for the disinfection of non-invasive medical device surfaces within medical and dental environments. Ideal to use where time is critical and in high traffic areas.

- Contact time of only 60 seconds
- Wide effective spectrum
- Virus tested
- Kills TB and Coronavirus
- Quick drying
- Leaves no smears or residue
- Dual registration; Medical Device and Biocide
- Fresh lemon aroma

Bossklein SDactiv is ideal to use on alcohol resistant, non-porous surfaces within the dental surgery. The 60 second contact time makes them perfectly suited where short contact times are important, eg in between patients.

The 63% alcohol based formulation has a powerful microbial action, the perfect solution for controlling infections on resistant surfaces in clinical environments.

SDactiv has been rigorously tested and proven to be effective against all Bacteria and Yeast. This highly effective formula has been tested to real world EN standards for proven effectiveness including Phase II Step I and II. Effective against Coronavirus, HIV, Hepatitis, Influenza, TB, Adeno and Noro viruses in just 60 seconds.

The fresh lemon aroma leaves equipment and rooms reassuringly fresh, clean and ready for the next patient.

COMPOSITION

Ethanol (64-17-5) 63% w/w (630g/kg), Aqua, Perfumes (Contains Linalool, Citral, Geraniol, Alpha-Isomethyl ionone).

INSTRUCTIONS FOR USE

Spray clean surfaces and objects to be disinfected from a distance of 30cm so that they are completely moistened and let the disinfectant take action for 60 seconds. Allow surface spray to evaporate. For professional use only.

REGISTRATION

Dual registered as a Class IIa medical device for use on both non invasive medical devices and as a general purpose Biocide registered disinfectant. This ensures suitability for use within all types of clinical and non clinical environments.

TEST STANDARDS

The following test standards have been achieved: EN14885, EN14476, EN14348, EN13727, EN13697, EN13624.

VARIANTS AND PACKAGING

BOS9782 500ml Trigger Spray x 10 per box
BOS9801 1 litre Trigger Spray x 10 per box
BOS9779 5 litre Refill Bottle x 2 per box

MICROBIOLOGICAL EFFICACY

MICROORGANISM	TIME
Bactericidal	60 sec
BS EN 13727:2012 + A2:2015	
EN 13697:2015+A1:2019 (E)	
Pseudomonas aeruginosa, Staphylococcus aureus, Enterococcus hirae, E. coli, VRE	60 sec
Yeast	60 sec
BS EN 13624:2013	
EN 13697:2015+A1:2019 (E)	
Candida albicans	60 sec
Enveloped viruses	60 sec
BS EN 14476:2013 + A1 2015	
Hepatitis B virus (HBV), Hepatitis C virus (HCV), HIV, Measles virus, Hepatitis Delta virus (HDV), Influenza virus, Coronavirus, Herpesviridae, Filoviridae (Ebola), Rubella virus, Rabies virus, Poxviridae	60 sec
Non-Enveloped viruses	60 sec
BS EN 14476:2013 + A1 2015	
Adenovirus 5, Murine norovirus	60 sec
Mycobactericidal	60 sec
BS EN 14348:2005	
Mycobacterium avium	60 sec
Tuberculocidal	60 sec
BS EN 14348:2005	
Mycobacterium terrae	60 sec

*Always test an inconspicuous area first.
Do not use on copper or aluminium.



Complete **hygiene** solutions

DirectaDentalGroup
DIRECTA TopDental ORSINO parkell KOHLER RÖNIGG polydentix
TrollDental Sendoline physics CONTACTZ ceramir

Topdental (Products) Ltd
Unit 1, Holmfield Industrial Estate
Halifax, West Yorkshire HX2 9TN
Tel +44 (0)1535 652 750 www.topdental.co.uk