FL-BOND II

en

Instructions for use

fr

Mode d'emploi

es

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ru

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LIGHT-CURE FLUORIDE RELEASING DENTAL BONDING SYSTEM

FL-BOND II forms part of a new category "GIOMER". Its one-bottle ethanol-based PRIMER is a further development of the previous 2-bottle PRIMER. The BONDING AGENT contains S-PRG filler, enhanced and based on the PRG technology. It provides a stable fluoride release and recharge and offers radiopacity.

PRECAUTIONS FOR DENTAL PERSONNEL AND PATIENTS

- Do not use this material on patients with known allergies or hypersensitivities to these types of chemicals. (This pertains to operators also.)
- If any eruption or skin inflammation occurs on a patient when using this product, immediately discontinue use and have the patient seek medical advice. (The same precaution is applied for the operators.)
- Avoid contacting soft tissues, skin or eyes. In case of contact, rinse immediately with copious amounts of water and if necessary, seek medical advice.

APPLICATION

- Bonding of composite resin to tooth structure
- Pre-treatment of resin-cement to tooth structure
- Repair of restorations

DIRECTIONS FOR USE

Bonding of composite resin (application 1) and pre-treatment of resin-cement to tooth structure (application 2)

- After cleaning tooth surface, take shade while the tooth is still moist, make cavity preparation, isolate prepared tooth with a rubber dam and protect the dentin close to the pulp with calcium hydroxide.
- Optional ETCHANT treatment

If the cavity is mostly enamel, the optional ETCHANT treatment would be preferable before priming. Apply the ETCHANT onto the whole cavity for 15 seconds, then rinse with water and air dry.
- Application of FL-BOND II

Apply PRIMER onto the whole cavity with MicroBrush and leave it for 10 seconds. Then, dry thoroughly with oil-free air. Apply BONDING AGENT evenly onto with MicroBrush. Do not blow air after applying BONDING AGENT. (Do not mix PRIMER ad BONGING AGENT in the dish.) If the layer of BONDING AGENT becomes too thin, the bonding strength decreases. Immediately light-cure the BONDING AGENT using light cure unit for 10 seconds (halogen lamp) or 5 seconds (high power type LED light).

4-1. [application 1] Complete the filling, light curing, and polishing, according to the instructions of the composite resin manufacturer.

4-2. [application 2] Complete adhesion according to the instructions of the resin-cement manufacturer.

Repair of restoration materials such as composite resin and crown and bridge resin materials (application 3)

- Pre-treat the surface to be repaired by optimizing the conditions depending on the material.
- Dispense BONDING AGENT onto the dish (do not mix BONDING AGENT with other liquids).
- Apply generously and evenly onto the whole cavity with McroBrush. Do not air blow.
- Immediately light-cure by using either a halogen lamp unit for 10 seconds or a LED light unit (high power type) for 5 seconds.
- Repair with composite resin. (Please follow the instructions of the manufacturer in the case of composite resin use.)

NOTE IN USE

- To avoid cross-infection, dispensed liquid should be used on one patient only and use chemical sterilization of the dish after each patient. Discard the use MicroBrush.

- Avoid using PREMIER and BONDING AGENT with other bonding materials.
- Wear of protective glasses when light curing is recommended of both operators and patients.
- The bottle should be tightly closed immediately after use.
- In case of bad feeling by inhaling vapor, immediately move to fresh air (both operators and patients).
- PRIMER and BONGIND AGENT are flammable. Keep away from source of ignition.
- If PRIMER of FL-BOND II contacts with oral soft tissue, the tissue may turn whitish, however, it will disappear in 2 or 3 days. (Don't touch the area.)
- Don't use this product of any other purpose except the purpose mentioned in the application and this product is intended for dental operators use only.
- Store this product in a refrigerator (1-10 °C/34-50 °F). Do not freeze. Bring it to room temperature before use. Use this product within the expiration date.

COMPOSITION

ETCHANT:	Water, Phosphoric acid
PRIMER:	Water, Ethanol, Carboxylic acid monomer, Phosphoric acid monomer and Initiator
BONDING AGENT:	S-PRG filler based on fluoro-boroaluminosilicate glass, UDMA, TEGDMA, 2-HEMA, Initiator

CAUTION: US Federal law restricts this device to sale by or on the order of a dental professional.

Location of the symbols glossary: https://www.shofu.com/en/symbols_glossary/

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