


NEW TePe EasyFit™



One brush for small and medium gaps

Because gaps between teeth vary, TePe EasyFit™ is a new tapered interdental brush that's designed to fit both small and medium gaps between teeth, making it easier for your patients to clean between teeth with just one brush.



SIMPLIFY CHOICE - ONE BRUSH FOR SMALL & MEDIUM GAPS

EasyFit™			S/M*						
ISO size (0-8)	0	1	2	3	4	5	6	7	8

*Indicative sizing vs ISO sizes for interdental brushes.

How to use TePe EasyFit™



Carefully insert the brush between your teeth until you feel resistance, and move back and forth.

For use between back teeth, curve the wire slightly.

Also suitable for braces.



Available now from HealthCare Essentials



178414 EasyFit 25pcs

178014 EasyFit 6pcs

TePe EasyFit™ S/M is suitable for small and medium gaps, comparable to Interdental Brush ISO size 2-4.