

NeoPUTTY®

PREMIXED ZERO-WASTE BIOACTIVE BIO CERAMIC



Positive displacement
plunger tip for zero waste.



~5 hrs
setting time*

*Sets in vivo in the presence of moisture



NeoPUTTY's ease of use offers predictable, consistent, quality results in a dry or wet environment, whether filling a perforation, performing vital pulp therapy, or during apical surgery.

—**Sayed Attar, DDS, MS**
Endodontist

PERFECT CONSISTENCY, NO MIXING, NO WASTE

NeoPUTTY's uniform, low-tack consistency from beginning to end, with a 3-year shelf life and no dry out between uses, ensures that no product goes to waste.

UNMATCHED RADIOPACITY

Highest in pre-mixed bioceramic putty class at 8.1 mm Al equivalent.

IMMEDIATELY WASH-OUT RESISTANT

Stays where you want it so you can efficiently complete the procedure ensuring bioactive benefits.