

Dryz®

Gingival Retraction Paste

Create a DRY Working Environment Every Time

Dryz® is compatible with the most popular retraction techniques which include compression caps, retraction cord, electrosurgery and lasers.

- ✓ Optimized for tissue management when taking impressions, restoring teeth with subgingival decay, seating restorations, placing rubber dam clamps, and bleaching teeth
- ✓ Quickly stops gingival bleeding and crevicular fluid seepage
- ✓ Easy clean-up, leaving no residue behind
- ✓ Set Time: 2 Minutes
- ✓ Available in two easily detectable contrasting colors

Take Control of the Sulcus for Picture Perfect Digital Impressions

When your prep is flawless, you want your digital scanner to easily capture every detail. Parkell's Dryz soft-tissue management paste helps you do that. This easy-to-use material fits comfortably into the workflow of every digital dentist, providing that extra push for clearing the field of unwanted obstacles, so you can quickly obtain an accurate and complete digital impression every time.



Green

Michael Karter, DDS



Blu



SKU#	DESCRIPTION
S183	Unit Dose Capsules Kit – Green ; (30) unit dose capsules (0.24g/capsule)
S180	Syringes Kit – Green ; (7) 0.5 ml (0.85g) syringes, (15) applicator tips
S182	Dryz Syringes VALUE PACK – Green ; (25) 0.5 ml (0.85g) syringes, (55) applicator tips
S186	Intro Kit – Blu ; (30) unit dose capsules (0.24g/capsule), Dispensing Gun
S187	Unit Dose Capsules Kit – Blu ; (30) unit dose capsules (0.24g/capsule)
S190	Syringes Kit – Blu ; (7) 0.5 ml (0.85g) syringes, (15) applicator tips
S191	Syringes VALUE PACK – Blu ; (25) 0.5 ml (0.85g) syringes, (55) applicator tips

ACCESSORIES & REFILLS

S181	Dryz Syringe Luer Lock Applicator Tips (40)
S189	Unit Dose Capsule Dispensing Gun



For more information, visit Parkell.com or contact your authorized Parkell distributor.

