

GamaCide³

Dental & Medical, Equipment DISINFECTANT CLEANER

- 3 Minute broad spectrum kill
- Disinfectant + cleaner
- Multi-surface
- Bactericidal, Fungicidal, Virucidal, & Tuberculocidal



COVER ALL YOUR SURFACES

Dental • Medical • Janitorial • Sanitation
EMS Veterinary • Pet Care Facilities
Body Art Studio • Athlete Centres • Spas

DIN 02358786

Product		GamaCide ³
Broad Spectrum	Bactericidal	✓
	Virucidal	✓
	Fungicidal	✓
	Tuberculocidal	✓
Surface Compliance	Plastic	✓
	Smooth/Textured Vinyl	✓
	Copper/Brass	✓
	Stainless Steel	✓
	Painted Surfaces	✓
	Plexiglas	✓
	Glass	✓
	Brushed Aluminum	✓
Safety	Biodegradable	✓
	Refillable Pouch	✓
	Health Rating "0"	✓
	Flammability "0"	✓
	Reactivity Rating "0"	✓
	pH	7.6

Classes of Microorganisms* | GamaCide³ Disinfection

HARDEST
to kill



GamaCide³

EASIEST
to kill

BACTERIA WITH SPORES

(Examples: Bacillus subtilis, Clostridium tetani, C. difficile, C. botulinum)

PROTOZOA WITH CYSTS

(Examples: Giardia lamblia, Cryptosporidium parvum)

MYCOBACTERIA

(Examples: Mycobacterium tuberculosis M. avium-intracellulare, M. chelonae)

NON-ENVELOPED VIRUSES

(Examples: Coxsackieviruses, polioviruses, rhinoviruses, rotaviruses, Norwalk virus, hepatitis A virus)

FUNGI

(Examples: Candida species, Cryptococcus species, Aspergillus species, Dermatophytes)

VEGETATIVE BACTERIA

(Examples: Staphylococcus aureus, Salmonella typhi, Pseudomonas aeruginosa, coliforms)

ENVELOPED VIRUSES

(Examples: Herpes simplex, varicella-zoster virus, cytomegalovirus, Epstein-Barr virus, measles virus, mumps virus, rubella virus, influenza virus, respiratory syncytial virus, hepatitis B and C viruses, hantaviruses, and human immunodeficiency virus)

*Source: Infection Control Guidelines: Hand Washing, Cleaning, Disinfection and Sterilization in Health Care (CCDR 24S8)

Surface friendly. Mild, neutral pH formula

To order GamaCide³ or for more information, contact your District Sales Manager or Germiphene at 1.800.265.9931 | info@germiphene.com | www.germiphene.com

GERMIPHENE®

GERMIPHENE®

A Division of Young Specialties

GamaCide³

Effectively kills the following microorganisms at room temperature: Mycobacterium bovis BCG, Duck Hepatitis B Virus (DHBV), Hepatitis C Virus (HCV), Rhinovirus ATCC VR-286, Poliovirus 1 ATCC VR-58, Herpes Simplex II Virus (HSV 2) G Strain, Human Immunodeficiency Virus (HIV-1), Influenza Virus Strain A2/Hong Kong, Staphylococcus aureus ATCC 6538, Staphylococcus aureus (MRSA), Enterococcus faecium (VRE), Pseudomonas aeruginosa ATCC 15442, Salmonella enterica ATCC 10708, Trichophyton mentagrophytes ATCC 9533.

Use on both hard and soft surfaces including: plastic, ceramic, stainless steel, chrome, acrylic, leather, vinyl, glass, Plexiglas®, and baked-on surfaces.

Compatible with most impression materials and dentures.

Non-sticky and leaves no hazardous chemical residue when used as directed.

Non-flammable. No special shipping or storage is required.

Safely dispose via standard waste disposal methods.

ACTIVE INGREDIENTS:

- Quaternary Ammonium Chloride.....**0.35% w/w**
- Isopropanol Alcohol.....**20.00% w/w**

REFILL
eco friendly
choice

Try both **liquid + wipes** to ensure all your infection control needs are met



10L



4L



160 regular wipes (15 x 17 cm)
& 160 refill pouch

To order GamaCide³ or for more information, contact your District Sales Manager or Germiphene at 1.800.265.9931 | info@germiphene.com | www.germiphene.com

GERMIPHENE®