



Dental Waterline Treatment Maintenance & Shock

Specifically designed for use in independent water bottle systems.

- Biofilm and bacteria control
- No harsh chemicals
- Odourless, colourless, tasteless, pH neutral
- Non-foaming, consistent results every time
- Tested on various waterline tubing with no negative effects or corrosion
- One E2O™ bottle lasts approximately 4–6 months per dental chair



Product Code:
GP-HE20

TAKE CHARGE of your office water today.
Easy dose pump. No measuring, no mess.

- Eco-friendly solution, E2O™ breaks down into water and oxygen
- Non-toxic to equipment or patients

- No waiting for tablets to dissolve, no residue/clogging from undissolved tablets
- No patient complaints about bad taste or grit from tablet residue




REGULAR TREATMENT



SHOCK TREATMENT



Water sample taken from an ultrasonic scaler



Before E2O™



After E2O™

Tests performed showed a significant reduction of bacteria in water sources after the use of E2O™.