

Tips on Tips

Universal Tips

Tissue Management Tips

Endodontic Tips

Periodontal Tips

Your Total Tip Solutions

Dental syringes are some of the most common pieces of equipment used by dental professionals in a wide range of dental procedures, from the most routine to the most complicated surgical procedures. At Vista Apex, we offer complete syringes, dental syringe tips (including bendable tips, dripless tips, brush tips, and more), and accessories like tip covers.

With many of these products, we offer the option to buy in bulk so your dental practice is always stocked with dental syringes and syringe tips. As with all our dental products, we aim to deliver the highest quality materials to make sure your practice is treating patients right.

SCAN HERE
For More Information



VISTA APEX
vistaapex.com
877.418.4782

Tissue Management

Stat-Flo™ Brush Tip

- Padded-end brush tip



20 cup
V-312102
100 cup
V-312100

Designed for:
Hemostatic agents

MST™ Universal Tip

- Material saving technology®



100 cup
22ga
V-315310

Designed for:
Flowable composites, pulp capping, etchants, & sealants

Spira-Flo™ Brush Tip

- Nylon brush fibers for controlled delivery
- Internal spiral
- Chamber minimizes bubbles



100 cup
V-321400

Designed for:
Etchants and resins

Sol-u-Flo™ Applicator Tip

- 19ga tip with adjustable brush fiber length
- Curved cannula for optimal entry



20 cup
V-312202
100 cup
V-312200

Designed for:
Aqueous solutions such as Caries Indicator, Silane, and Chlorhexidine

Exacta-Flo™ Brush Tip

- Adjustable, soft brush fibers minimize bubbles



100 cup
V-321610

Micro-Flo™ Tip

- Polyimide tip for microscope enhanced dentistry



50 cup
26ga
V-406626
27.5ga
V-406627

Designed for:
Flowable composites, pulp capping, etchants, & sealants

Volu-Flo™ Tip

- Delivers large volumes of material (18ga)



100 cup
V-321100

Designed for:
Viscous materials, pastes & putties

1/2" Bendable Applicator Tip

- Blunt end, dead-soft needle resists crimping



22ga
20 cup
V-315202
100 cup
V-315222



25ga
20 cup
V-315205
100 cup
V-315225

Designed for:
Etchants, resins, and flowable composites

Blue-Flo™ Tip

- 23ga soft plastic tip



20 cup
V-316202
100 bag
V-316222

Designed for:
Irrigation, aeration & impression materials

Appli-Vac™ 3/4" Bendable Tip

- Blunt end, dead-soft needle resists crimping



100 cup
16ga
V-315116



100 cup
18ga
V-315118



100 cup
19ga
V-315119



100 cup
20ga
V-315120



100 cup
22ga
V-315122



100 cup
25ga
V-315125

Designed for:
Sealants, flowable composites, etchants, resins, primers

Pre-Bent Tip

- Blunt end pre-bent needle
- Opaque photo sensitive hub (22/25ga)



100 cup
16ga
V-315316



100 cup
20ga
V-315320



100 cup
22ga
V-315322



100 cup
25ga
V-315325

Designed for:
Etchants, resins and flowable composite

MST Visco-Tip™

- Bendable, flexible 25ga tip




20 cup
V-312302
50 cup
V-312300

Designed for:

Viscous materials such as bioceramic sealers and calcium hydroxide

Voyager Irrigating Tip

- Dual side ports for safety and efficacy
- Enhanced flexibility to navigate canals
- Silicone stopper for depth reference



17mm 27ga V-315727	17mm 29ga V-315729	17mm 30ga V-315730	17mm 31ga V-315731
27mm 27ga V-315827	27mm 29ga V-315829	27mm 30ga V-315830	27mm 31ga V-315831

Designed for:

Endo irrigating solutions

Appli-Vac™ 1" Bendable Irrigating Tips

- Slotted and side-vented




20 cup 23ga V-315023	20 cup 27ga V-315027	20 cup 30ga V-315030
100 bag 23ga V-315523	100 bag 27ga V-315527	100 bag 30ga V-315530

Designed for:

Endo Irrigating solutions (27 ga / 30 ga)
Perio Irrigation (23 ga)

Secure-Lock™ Appli-Vac™ 1" Irrigating Tip

- Slotted and side-vented
- Secure-Lock™ threads lock tip in place



80 bag 23ga V-314523	80 bag 27ga V-314527	80 bag 30ga V-314530
-----------------------------------	-----------------------------------	-----------------------------------

Designed for:

Endo Irrigating solutions (27 ga / 30 ga)
Perio Irrigation (23 ga)

Vista-Probe™ Irrigating Tip

- Closed-end for side-port delivery
- Lateral delivery of irrigants reduces apical pressure



20 cup 23ga V-315123	20 cup 27ga V-315127	20 cup 30ga V-315130
100 bag 23ga V-315623	100 bag 27ga V-315627	100 bag 30ga V-315630

Designed for:

Endo Irrigating solutions (27 ga / 30 ga)
Perio Irrigation (23 ga)

Secure Lock™ Vista-Probe™ Tip

- Closed-end for side-port delivery
- Lateral delivery of irrigants reduces apical pressure
- Secure-Lock™ threads lock tip in place



80 bag 23ga V-314623	80 bag 27ga V-314627	80 bag 30ga V-314630
-----------------------------------	-----------------------------------	-----------------------------------

Designed for:

Endo Irrigating solutions (27 ga / 30 ga)
Perio Irrigation (23 ga)

Elasti-Vac™ Tip

- Narrow, flexible tip easily navigates curved canals



50 cup
V-406150

Designed for:

Endodontic evacuation

Versa-Brush™

- Rotational utility brush



40 cup
V-204300

Designed for:

Cleaning post holes and coronal access openings

NiTi SuperFlex™ Tip

- AUTOCLAVABLE. Slotted and side-vented, 30ga for safe irrigation
- 2.5x more flexible



17mm (short) 6 cup V-405107	25mm (long) 6 cup V-405108	50 cup V-406350
Mixed 6 cup V-405106		

Designed for:

Endo irrigating solutions

Flexi-Glide™ Utility Tip

- 30ga flexible tip for easy access

Periodontal

Pocket-Probe™

- Closed-end with multiple perforations for thorough flushing



100 cup
21ga
V-316102

Designed for:

Flushing of perio pockets or inter-dental spaces