PARADISE
DENTAL TECHNOLOGIES

Introducing the new
Boge Precision Duo™
R117



a PDT signature anterior design combining a **modified hoe scaler**, built for ergonomics, control, and versatility, with a **Nebraska 128** thin sickle scaler for optimal supragingival clinical applications.



- Reinforced terminal shank design resists flexion for confident removal of calculus, stain, and residual cement.
- Designed for compatibility with the Gleason Guide™ and other sharpening systems for easy maintenance and resharpening
- Engineered with PDT metallurgy and a lightweight resin handle for lasting sharpness, superior grip, and exceptional clinical performance.

Created in collaboration with
Dr. Emily Boge ,EdD, RDH, CDA



(406)626-4153
customer.care@pdtdental.com
pdtdental.com